


Present


- Classification : Heart, Respiration sound
- Heart : Normal, Abnormal
- Lung : Normal, Crackle, Wheezing, Rhonchi, Stridor

Present


- Quality :
 - ✓ Sensitivity 99%
 - ✓ Specificity 99%
 - ✓ Accuracy 89%



System configuration of JABES-AI



AI Server



Use's smart phone




JABES-Android App

JABES-Apple App



Doctor's Clinic






System configuration of JABES-AI



+ JABES App for Android, Apple

The JABES-Handy connects to a smart phone

Preparation



After connecting the stethoscope and the enclosed cable to the smartphone, enter the JABES app in web and use it according to the manual.